

ENDO LINEAR CUTTER STAPLER WITH NON-STOP ARTICULATING LEVER



Extra Grip Reloads

- Atraumatic extra grip
Tissue affectively held in the jaws
- Reduced firing force
- Uniform staple line
- Less slippage



Reloading unit

Load alignment indicator:
Precise loading
of the reloading unit
into the handle



Curved tip Reloads
specially indicated for
thin tissue and vessels



The most complete range of endo reloads
to cover all surgery procedures



Five different staple heights,
2.5mm, 3.5mm, 4.2mm, 4.8mm, 5.0 mm,
addressed to various tissue thicknesses

Universal handle for all reloading units, both straight and articulating



Ergonomic
Non-slipping Handle



Innovative 45 degrees
Non-Stop articulation system

Universal handle with
better design and ergonomics
for easy 360 degrees knobbing

NEW GENERATION ENDO STAPLING SYSTEM



Product Name	Product	Description	Color
UNIVERSAL Endo Linear Cutter	CEAB-30	Endo Linear Cutter with Non-Stop articulation system (Standard Handle)	
	CEAB-45	Endo Linear Cutter with Non-Stop articulation system (Short Handle)	
	CEAB-60	Endo Linear Cutter with Non-Stop articulation system (Extra Long Handle)	

Product Name	Product	Description	Color
Endo Articulating Reloading	CADF-30T	PANTHER Endo Articulating reloading unit 30 mm with P2G Technology and 2.5mm staple height (White)	Articulated reload ○ White
	CADF-30N	PANTHER Endo Articulating reloading unit 30 mm with P2G Technology and 3.5mm staple height (Blue)	Articulated reload ● Blue
	CADF-45T	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 2.5mm staple height (White)	Articulated reload ○ White
	CADF-45N	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 3.5mm staple height (Blue)	Articulated reload ● Blue
	CADF-45S	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 4.2mm staple height (Gold)	Articulated reload ● Gold
	CADF-45D	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 4.8mm staple height (Green)	Articulated reload ● Green
	CADF-45R	PANTHER Endo Articulating reloading unit 45 mm with P2G Technology and 5.0mm staple height (Black)	Articulated reload ● Black
	CADF-60T	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 2.5mm staple height (White)	Articulated reload ○ White
	CADF-60N	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 3.5mm staple height (Blue)	Articulated reload ● Blue
	CADF-60S	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 4.2mm staple height (Gold)	Articulated reload ● Gold
	CADF-60D	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 4.8mm staple height (Green)	Articulated reload ● Green
	CADF-60R	PANTHER Endo Articulating reloading unit 60 mm with P2G Technology and 5.0mm staple height (Black)	Articulated reload ● Black

New Generation Endo Stapling System

P2G



Panther Healthcare Medical Equipment Co., Ltd
Beijing | New York | São Paulo | Toronto | Mumbai

Tel: +8610-80789058
www.pantherhealthcare.com
Email: service@pantherhealthcare.com

www.pantherhealthcare.com

Endo Linear Cutter Stapler

WITH NON-STOP ARTICULATING LEVER

