



**BEAUTY COUCH  
PODIATRY CHAIR  
FACIAL BED  
TREATMENT TABLE**



# COMPANY OVERVIEW



Tronwind Industries Co. Limited is a leading supplier of medical chairs who is located in Southern China. We provide a comprehensive range of, beauty treatment beds, chairs and stools which are essential for any doctors, dentists, nurses, technicians, patients and office workers to perform effectively and comfortably. Our main products are doctor stools, nurse stools, saddle stools, laboratory chairs, beauty beds, facial beds, podiatry chairs, treatment beds and office furniture with full range.

We try our best to stay on top of the industry, provide exceptional value and high quality products. Our ergonomically designed chairs, reduces back and shoulder fatigue, especially for doctors, patients and operators, who have to sit for long periods, or in specific positions. Our chairs are typically used in all therapy areas, such as Hospital, Clinic, Theatre, Maternity, Dentistry, Laboratories, Beauty Salon, Massage Center and many more clinical departments, throughout hospitals from all over the world.

We are seeking partnerships and distributors to promote & sell our range of chairs. Therefore, traders, distributors are welcome. Also we do OEM for customers. Tronwind is committed to giving you great product selection, good prices, and world class customer service!



# TRA03

## ELECTRIC BEAUTY CHAIR PODIATRY CHAIR / BEAUTY BED TREATMENT COUCH

With attractive design, Tronwind TRA03 Beauty Chair / Multifunction Treatment Bed allow independent adjustments that make treatment easier and more comfortable for both patients and staff. The wide range of leg rest adjustments provides a perfectly vertical upright sitting position. Fully height adjustable, offers generous seat tilt and backrest angle change all by remote control. The solid upholstered armrest makes it possible to achieve the most suitable and restful arm position for any physique.



45°+45°  
or 135°+45°

**4  
MOTORS**

CE



### LEATHER COLORS



### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment



Footrest-adjustment



Tilt-adjustment

### SPECIFICATIONS

Length	190cm	Width	60cm
Max. Height	92cm	Min. Height	61cm
Max. Load	250kg	Armrest	50cm
Net Weight	95kg	Gross Weight	105kg
Motors	4 or 3 Motors (CE)	Packing Size	178×84×77cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	Fiber Leather	Control	Handset Control/Foot Control (Optional)





# TRA02

## ELECTRIC PODIATRY CHAIR PEDICURE CHAIR / BEAUTY BED OPERATING COUCH

With attractive design, Tronwind TRA02 Podiatry Chair / Multifunction Treatment Bed allow independent adjustments that make treatment easier and more comfortable for both patients and staff. The wide range of leg rest adjustments provides a perfectly vertical upright sitting position. Fully height adjustable, offers generous seat tilt and backrest angle change all by remote control. The solid upholstered armrest makes it possible to achieve the most suitable and restful arm position for any physique.



45°+45°  
or 135°+45°

**5  
MOTORS**  
CE



### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment



Footrest-adjustment



Tilt-adjustment

### SPECIFICATIONS

Length	193cm	Width	60cm
Max. Height	92cm	Min. Height	61cm
Max. Load	250kg	Armrest	50cm
Net Weight	95kg	Gross Weight	105kg
Motors	5 or 3 Motors (CE)	Packing Size	178×84×77cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	Fiber Leather	Control	Handset Control/Foot Control (Optional)

### LEATHER COLORS





# TRW03


## TREATMENT CHAIR BEAUTY CHAIR / FACIAL BED HAIR TRANSPLANT CHAIR

Tronwind Electric Treatment Chair TRW03 is the ultimate in versatility and comfort. This is a hair transplant chair/derma chair with motorized positioning control that can be adjusted to different positions through a handset remote controller. It can be adjusted individually in terms of backrest, footrest seat tilt and height. Up-Down movements of the armrests are operated vertically. This electric treatment chair is designed ergonomically according to the needs of customer in hair transplantation, medical treatment and beauty treatment.


### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment



Height-adjustment



Tilt-adjustment

### SPECIFICATIONS

Length	195cm	Width	62cm
Max. Height	88cm	Min. Height	68cm
Max. Load	250kg	Armrest	40cm
Net Weight	90kg	Gross Weight	100kg
Motors	4 or 3 Motors (CE)	Packing Size	178×84×80cm
Upholstery	High Density Foam	Base	Stainless Steel
Leather	Fiber Leather	Control	Handset Control/Foot Control (Optional)



### LEATHER COLORS



4  
MOTORS

CE



# TRW01

## POWER PROCEDURE CHAIR BEAUTY CHAIR / FACIAL BED TREATMENT BED

Tronwind Power Procedure Chair TRW01 is designed for practically all kinds of beauty treatments/cosmetology treatments as well as other massage and spa treatments. Tronwind Power Procedure Chair TRW01 has ergonomic and stylish design and excellent technical features. Enables specialists to work with any part of patient's body easily. The angles of the bed have soft round shape and are absolutely safe. Comfortable, simple and user-friendly functions of this model provide patients with maximum comfort and therapeutic effect.



**4  
MOTORS**  
CE

### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment



Height-adjustment



Tilt-adjustment

### SPECIFICATIONS

Length	190cm	Width	60cm
Max. Height	84cm	Min. Height	64cm
Max. Load	250kg	Armrest	50cm
Net Weight	80kg	Gross Weight	90kg
Motors	4 or 3 Motors (CE)	Packing Size	176×75×74cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	Fiber Leather	Control	Handset Control/Foot Control (Optional)

### LEATHER COLORS





# TRW02

## TREATMENT CHAIR BEAUTY CHAIR / FACIAL BED MEDICAL COUCH

Tronwind Electric Treatment Chair TRW02 is designed for practically all kinds of beauty treatments/cosmetology treatments as well as other massage and spa treatments. Tronwind Treatment Chair TRW02 has ergonomic and stylish design and excellent technical features. Enables specialists to work with any part of patient's body easily. The angles of the bed have soft round shape and are absolutely safe. Comfortable, simple and user-friendly functions of this model provide patients with maximum comfort and therapeutic effect.



4  
MOTORS

CE

### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment



Height-adjustment



Tilt-adjustment

### SPECIFICATIONS

Length	190cm	Width	60cm
Max. Height	85cm	Min. Height	60cm
Max. Load	250kg	Armrest	50cm
Net Weight	80kg	Gross Weight	90kg
Motors	4 or 3 Motors (CE)	Packing Size	160×70×76cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	Fiber Leather	Control	Handset Control/Foot Control (Optional)

### LEATHER COLORS





# TRW04

## GYNECOLOGICAL EXAMINATION CHAIR GYNECOLOGY TABLE / MEDICAL BED TREATMENT CHAIR

TRONWIND gynecological examination chair with electric height adjustment by one actuator motor, includes gynecological leg holders, removable tray and ergonomic cushion. This gynecology table with 3 sections, electrically adjustable height controlled by a control panel. Backrest and legrest with pneumatic adjustment by gas strut. The upholstery are made with fire-retardant flexible soft high density foam filling covered with the washable fire-retardant PU leather, it is ergonomic shaped and there are fully adjustable backrest, legrest, leg holders.

### ADJUSTMENTS



Backrest-adjustment



Legrest-adjustment




Height-adjustment



Leg Supports Height adjustment

### SPECIFICATIONS

Length	186cm	Width	60cm
Max. Height	72cm	Min. Height	52cm
Max. Load	250kg	Armrest	50cm
Net Weight	85kg	Gross Weight	95kg
Motors	1 Motor (CE)	Packing Size	160×75×82cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	PVC Leather	Control	Handset Control/Foot Control (Optional)



### LEATHER COLORS



1  
MOTOR

CE



# TEP05

## ELECTRIC BEAUTY CHAIR BEAUTY BED / FACIAL BED MASSAGE TREATMENT BED

Tronwind electric beauty bed TEP05 is designed for practically all kinds of beauty treatments/cosmetology treatments as well as other massage and spa treatments. Tronwind electric beauty bed TEP05 has ergonomic and stylish design and excellent technical features. Enables specialists to work with any part of patient's body easily. The angles of the bed have soft round shape and are absolutely safe. Comfortable, simple and user-friendly functions of this model provide patients with maximum comfort and therapeutic effect.



**4  
MOTORS**

CE

### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment



Height-adjustment



Tilt-adjustment

### SPECIFICATIONS

Length	190cm	Width	62cm
Max. Height	85cm	Min. Height	60cm
Max. Load	250kg	Armrest	53cm
Net Weight	76kg	Gross Weight	85kg
Motors	3 or 4 Motors (CE)	Packing Size	145×66×76cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	Fiber Leather	Control	Handset Control/Foot Control(Optional)

### LEATHER COLORS





# TEP01

## ELECTRIC BEAUTY SPA BED BEAUTY COUCH MASSAGE TREATMENT BED

Tronwind motorized beauty bed offers soft and supportive PU upholstery to ensure your clients stay comfortable and relaxed, even for longer time treatments. The beauty bed has 4 motors to control the backrest, the height, the leg rest and the inclination. Both the head section and the foot section are extensible, or alternatively can be removed. The vent hole is hidden under the hole cushion of the round face and the arms can be removed or flipped back to allow you closer access to your customer or use it for comfortable arms.

### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment



Height-adjustment



Tilt-adjustment

### SPECIFICATIONS

Length	190cm	Width	62cm
Max. Height	92cm	Min. Height	68cm
Max. Load	250kg	Armrest	50cm
Net Weight	76kg	Gross Weight	85kg
Motors	3 or 4 Motors (CE)	Packing Size	145×66×76cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	Fiber Leather	Control	Handset Control/Foot Control (Optional)



### LEATHER COLORS



4  
MOTORS  
CE





# TEP02

## ELECTRIC PODIATRY CHAIR PEDICURE CHAIR / BEAUTY BED OPERATING COUCH

With attractive design, Tronwind TEP02 Podiatry Chair / Multifunction Treatment Bed allow independent adjustments that make treatment easier and more comfortable for both patients and staff. The wide range of leg rest adjustments provides a perfectly vertical upright sitting position. Fully height adjustable, offers generous seat tilt and backrest angle change all by remote control. The solid upholstered armrest makes it possible to achieve the most suitable and restful arm position for any physique.

### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment



Height-adjustment



Tilt-adjustment

5  
MOTORS

CE



### SPECIFICATIONS

Length	190cm	Width	62cm
Max. Height	92cm	Min. Height	68cm
Max. Load	250kg	Armrest	50cm (Wooden Cover Optional)
Net Weight	79kg	Gross Weight	88kg
Motors	4 or 5 Motors (CE)	Packing Size	145×66×76cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	Fiber Leather	Control	Handset Control/Foot Control(Optional)

### LEATHER COLORS





# TEP03

## MOTORIZED FACIAL BED BEAUTY COUCH / MASSAGE TABLE COSMETIC CHAIR

Tronwind motorized facial bed TEP03 offers special support backrest that provide extra lumbar support and good fit to ensure your clients stay comfortably during treatment with the easy-to-clean PU upholstery. The 4-motors design allows you to control the backrest, the height, the leg rest and the angle of bed. The tapered headrest can be rotated to opposite or removed, which allows you can get closer to your client to do facial treatment. In addition, the inclinable facial bed is ideal for facial treatment or for injectable treatment.

### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment



Height-adjustment



Tilt-adjustment

### SPECIFICATIONS

Length	195cm	Width	62cm
Max. Height	83cm	Min. Height	60cm
Max. Load	250kg	Armrest	50cm (Wooden Cover Optional)
Net Weight	78kg	Gross Weight	88kg
Motors	3 or 4 Motors (CE)	Packing Size	166×67×69cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	Fiber Leather	Control	Handset Control/Foot Control(Optional)



### LEATHER COLORS



4

MOTORS

CE



# TEP04

## TREATMENT BED MASSAGE BED / PODIATRY CHAIR BEAUTY CHAIR / EXAM TABLE

Tronwind Company is proud to announce the availability of the Euro Style Beauty Massage Bed / Multi-Function Treatment Bed TEP04. This state of the art motorized chair offers clients unparalleled comfort and nail technicians' ergonomic, time-saving versatility.

### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment



Height-adjustment



Tilt-adjustment

### SPECIFICATIONS

Length	202cm	Width	62cm
Max. Height	92cm	Min. Height	68cm
Max. Load	250kg	Armrest	50cm (Wooden Cover Optional)
Net Weight	82kg	Gross Weight	92kg
Motors	3 or 4 Motors (CE)	Packing Size	168×69×76cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	Fiber Leather	Control	Handset Control/Foot Control(Optional)



### LEATHER COLORS



4  
MOTORS

CE





# TEP06

## ELECTRIC BEAUTY COUCH FACIAL BED / TREATMENT COUCH MASSAGE TABLE

The Tronwind electric beauty couch can do it all! Ideal for all beauty treatments, facial treatment, massage treatment, Tronwind treatment table can go from a chair – great for face treatments and makeup to a flat table suitable for massage, facials and advanced beauty treatments. The 4 motors make the table able to adjust in all directions including upright, flat and tilt. This Tronwind Beauty Couch TEP06 allows your clients to be comfortable while receiving the ultimate indulgence or relaxation treatment.

### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment



Height-adjustment



Tilt-adjustment

### SPECIFICATIONS

Length	190cm	Width	62cm
Max. Height	83cm	Min. Height	60cm
Max. Load	250kg	Armrest	50cm
Net Weight	78kg	Gross Weight	88kg
Motors	3 or 4 Motors (CE)	Packing Size	150×67×69cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	Fiber Leather	Control	Handset Control/Foot Control(Optional)



### LEATHER COLORS



4  
MOTORS

CE



# TAE05

## MULTI PURPOSE MASSAGE TABLE ELECTRIC SPA BED / THERAPY BED ESTHETICIAN BED / MEDICAL COUCH

Tronwind TAE05 Electric Multi Purpose Massage Table with 4 motors for body, face and massage treatments. This massage treatment bed provides maximum comfort for both client and practitioner, thanks to 4 actuator motors for height, backrest and legrest adjustments. The durable and soft polyurethane upholstery is easy to keep clean with a warm cloth and soapy water or leather wax. The adjustable armrests and breath hole provide ultimate client comfort. Two pillars supports and strong steel structure guarantee stability and safety of this massage table when heavy loading.



**4  
MOTORS**  
CE


### ADJUSTMENTS



Armrest-adjustment



Legrest-adjustment



Height-adjustment



Backrest-adjustment

### SPECIFICATIONS

Length	203cm	Width	76cm
Max. Height	83cm	Min. Height	60cm
Max. Load	250kg	Armrest	Adjustable
Net Weight	100kg	Gross Weight	115kg
Motors	4 Motors (CE)	Packing Size	205×80×80cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	PVC Leather	Control	Handset Control/Foot Control (Optional)

### LEATHER COLORS





# TAE04

## ELECTRIC SPA BED ELECTRIC MASSAGE BED / THERAPY BED ESTHETICIAN BED / MEDICAL COUCH

Tronwind Electric Spa Bed TAE04 is designed for ultimate comfort. It's comfortable, reliable, easy to use and with the combination of elegance and simplicity. This is 3-motor SPA bed that adjust the height, the backrest and legrest (center). 2-motor and 1-motor versions are optional. With its three quiet actuator motors, ensures relaxation, every adjustment is smooth and quiet. PU upholstery with high quality and comfort padding. Headrest with breather hole.



### LEATHER COLORS



### ADJUSTMENTS



Backrest-adjustment



Height-adjustment



Legrest-adjustment

### SPECIFICATIONS

Length	200cm	Width	67cm
Max. Height	84cm	Min. Height	59cm
Max. Load	250kg	Base	120×58cm
Net Weight	100kg	Gross Weight	110kg
Motors	3 Motors (CE)	Packing Size	205×77×80cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	PVC Leather	Control	Handset Control/Foot Control (Optional)





# TAE01

## ADJUSTABLE MASSAGE BED TREATMENT TABLE / THERAPY BED OPERATING BED / MEDICAL COUCH

The Tronwind Massage Bed / Beauty Bed TAE01 features an electric height adjustment conveniently operated by a hand control. The quiet motor effortlessly raises and lowers the table. The solid steel frame ensures stability with the added bonus of electric back and leg rest adjustment, which can be easily converted into a chair position. With the removable flexible armrests, Tronwind Massage Bed / Beauty Bed TAE01 offers incomparable quality and value for beauty treatments, facial treatment, spa, and massage and clinic practitioners.




**2  
MOTORS**  
CE

### ADJUSTMENTS



Backrest-adjustment



Height-adjustment

### SPECIFICATIONS

Length	190cm	Width	67cm
Max. Height	79cm	Min. Height	60cm
Max. Load	250kg	Base Size	114×56cm
Net Weight	95kg	Gross Weight	110kg
Motors	2 Motors (CE)	Packing Size	195×72×68cm
Upholstery	High Density Foam	Base	Stainless Steel
Leather	Fiber Leather	Control	Handset Control/Foot Control (Optional)

### LEATHER COLORS





# TEP07

## COSMETIC BED / BEAUTY CHAIR FACIAL BED / SPA TREATMENT BED


The Tronwind Cosmetic Bed/Beauty Bed TEP07 features an electric height adjustment conveniently operated by a hand control. The quiet motor effortlessly raises and lowers the table. The solid steel frame ensures stability with the added bonus of electric back and leg rest adjustment, which can be easily converted into a chair position. With the removable flexible armrests, Tronwind Cosmetic Bed/Beauty Bed TEP07 offers incomparable quality and value for beauty treatments, facial treatment, spa, and massage and clinic practitioners.



### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment



Height-adjustment

### SPECIFICATIONS

Length	190cm	Width	62cm
Max. Height	80cm	Min. Height	58cm
Max. Load	250kg	Armrest	72cm
Net Weight	76kg	Gross Weight	85kg
Motors	3 Motors (CE)	Packing Size	160×67×70cm
Upholstery	High Density Foam	Base	Steel Structure, Reinforced Plastic Casing
Leather	Fiber Leather	Control	Handset Control/Foot Control(Optional)

### LEATHER COLORS





# TEP08

## ADJUSTABLE BEAUTY TREATMENT CHAIR TREATMENT TABLE / THERAPY BED OPERATING BED

Tronwind Electric Treatment Bed / Beauty Massage Bed TEP08 have a very unique design with unlimited positioning combination which will give you the ultimate freedom to perform even outpatient aesthetic procedures and other treatments. This is ultimate luxury facial and massage bed with ultra-luxurious breathable fabric which is very easy to clean and maintain. Stylish hand held remote with long cord. Even at the steepest angel the bed will maintain its excellent weight distribution system. With the option of removable flexible armrests,


### ADJUSTMENTS



Backrest-adjustment



Footrest-adjustment




Height-adjustment



Tilt-adjustment

### SPECIFICATIONS

Length	190cm	Width	62cm
Max. Height	88cm	Min. Height	68cm
Max. Load	250kg	Armrest	50cm
Net Weight	82kg	Gross Weight	92kg
Motors	3 or 4 Motors (CE)	Packing Size	160×67×76cm
Upholstery	High Density Foam	Base	Stainless Steel 85×53cm
Leather	Fiber Leather	Control	Handset Control/Foot Control(Optional)



### LEATHER COLORS





# TB01 / TB02 / TB03

THERAPIST STOOL / SADDLE STOOL  
SALON STOOL / DOCTOR STOOL



**TB 01**



**TB 03**



**TB 02**

## ADJUSTMENTS



Backrest tilt adjustment



Seat height adjustment



Backrest height adjustment

## SPECIFICATIONS

Backrest	L340mm×W410mm	Seat	L420mm×W380mm
Gas Lift	140mm travel distance	Base	Φ560mm Aluminium alloy
Castor	Φ50mm	N.W.	9 kg
Upholstery	High density foam	Cover	Fiber leather

## LEATHER COLORS







Barosupport.com  
@barosupportrd

info@barosupport.com.do

809-732-4440

📞 809-769-6025

